PRODUCED BY:
KONUS ITALIA GROUP s.r.l.
 Via Mirandola, 45 - 37026 SETTIMO di P. (Verona)
 ITALY - Tel. +39-45-6767670 • Fax +39-45-6767671

KONUS IS A TRADE MARK AND REGISTERED
 OF KONUS ITALIA GROUP s.r.l. - ITALY
 & SELECTED COUNTRIES

Branch offices:
KONUS FRANCE s.a.l.
 80, Rue Charles Duschesne - "Le Mercure C"
 13851 AIX EN PROVENCE Cedex 3 - FRANCE Tel.
 +33-442-39 41 78 • Fax +33-442-39 77 12

KONUS ESPAÑA S.L.
 C/Caucho, 53-28850 - TORREJON DE ARDOZ
 (Madrid) - ESPAÑA
 Tel. +34-91-67 78 663 • Fax +34-91-67 78 538

BIOLOGIC BINOCULAR MICROSCOPE KONUS ACADEMY

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ENGLISH

We appreciate your choice on our Konus Microscope.

BIOLOGIC BINOCULAR MICROSCOPE KONUS ACADEMY

TECHNICAL CHARACTERISTICS

- Biological monocular microscope with pivoting head a 360°
- Double focusing device (macrometric and micrometric)
- N.A. 1.2 condenser with iris diaphragm and filterholder.
- 10x eyepiece
- Illumination with halogen lamp 12V/10W (2 pieces) and external transformer.
- Four achromatic objectives (4x; 10x; 40x retractable and 100x immersion)
- Magnifications in equipment: 40x; 100x; 400x; 1000 times.
- Movable stage.
- Coloured filters.
- Anti-dust cover.
- 2 years warranty.

COMPONENTS (see figure)

1. Eyepiece
2. Head
3. Stand
4. Base
5. Rotating turret objective holder.
6. Objectives.
7. Movable stage.
8. Specimen table.
9. Condenser with iris diaphragm and filter holder.
10. Halogen lamp (internal)
11. On/off switch.
12. Macrometric focusing.
13. Micrometric focusing (fine regulation).

MOUNTING.

- Open the packaging and take out all the components.
- Position the base (4) with the stand (3) and head (2) on a dry and stable table.
- Take the head (2) and insert it on the upper part of the stand (3).
 - Take the 10x eyepiece (1) and after taking off the caps, insert them into the head (2).
 - Take out the Objectives (6) from their casing and screw them onto the turret (5), turn the turret until the 4x objective (the shortest) is directly above the specimen stage (8).
 - Take the transformer and plug it into an electric socket, make sure that the data on the label of the transformer matches that of your electrical supply).

- Insert the jack plug of the transformer into the socket of the base (4) of the instrument.
- Turn on the light with the two way switch (11). Now your microscope is ready for use.

CHANGING THE BULB: The Konus Academy model is supplied with a Halogen Lamp 12V/10W. In case it needs replacing, take the following steps:

1. Remove the jack plug from the microscope.
2. Turn the micrometric focus knob (13) lifting the specimen table (8) until the limit.
3. Move the condenser (9) as high as possible.
4. Unscrew the halogen lamp holder.
5. Replace the lamp with a fresh one.
6. Screw on back the lamp

40x to 100x are much more longer than the objectives of 4x and 10x, and so, rotating the objective turret can touch the specimen slide (8), in this case, it is necessary to first lower the specimen stage, so as not to create any damage.

You will notice that during focusing of the 40x and 100x objectives, they will first touch the cover slide and then the terminal section will enter the inside of the objective. This is absolutely necessary for high enlargements and they have been prepared for such operations with a special spring.

Be very careful not to put too much pressure onto the focusing knob as the specimen slide might be cracked by the force.

The 100X objective supplies a high number of enlargements, and in these conditions, the air present around the objective can highly worsen the quality of the image. To take care of this problem, an immersion oil may be used (optional, art. 5182); using one or two drops of the 5182 oil on top of the cover slide, so as not to have any air between the objective and the specimen during observation.

During the observation of very transparent samples, the light might be too intense to be able to see all the details. To resolve this problem, the condenser must be adjusted (9) with an iris diaphragm or coloured filters.

Highering up the condenser and/or opening the diaphragm means letting in more light, while lowering the eyepiece and/or closing the diaphragm means decreasing the light. For particular specimens it is possible to use appropriate dyes to bring out the details.

MAINTENANCE: The instrument does not present itself with particular maintenance problems, it is important to follow these simple rules:

Do not tamper with the instrument, especially the electrical wiring and the optical parts. Protect the instrument from dust, acids, shocks and/or impacts, water and/or humidity. Do not place the transformer within reach of children.

The eventual vision by children of the microscope should be taken under the supervision of adults.

OPTIONAL ACCESSORIES: Photo adapter #5126 (1/2 ring not included) - eyepiece 15x #5115 - CCD camera adapter #1599 - CCD camera #5821 - replacement bulb #5234.

PREPARED SLIDES: Code **4855** Fish Frogs And Other Amphibians. - Code **4856** Lizards, Snakes And Birds. - Code **4862** Animal Reproduction. - Code **4863** Animal Embryology And Development. - Code **4877** Bacteria, Simple Organisms. - Code **4880** Tropical Plants. - Code **4881** Plant Reproduction And Propagation. - Code **4886** The Wonderful Structure Of Vegetable Cells.